

# Usability Test Results

**Siemens ACUSON S3000™  
Ultrasound System, HELX™ Evolution  
with Touch Control**



## Research Goal

Siemens Medical Solutions USA, Inc. Ultrasound Business Area made a concerted effort to improve the usability of their family of premium ultrasound systems, the new ACUSON S Family™ of ultrasounds systems, HELX™ Evolution with Touch Control.

Siemens felt confident that they had achieved their goals and wanted Macadamian, a user experience design and development firm, to conduct an independent usability test to measure user performance and validate user satisfaction with one of these systems, the ACUSON S3000™ Ultrasound System.

# Macadamian Technologies (1)



Macadamian provides a complete range of high-quality usability, design and software engineering services to industry leaders across North America. Macadamian's work is founded in user-centered design to deliver context-aware and adaptive experiences.


For over 10 years, Macadamian's User Experience team has conducted usability assessments to better understand the experience users have or want to have with a particular product. These assessments include expert reviews, concept design walkthroughs, and usability testing on new and existing products for a wide variety of industries.

# Macadamian Technologies (2)



# Usability Test

This test is a user-centered design methodology that evaluates a product or service by testing it with representative users. Users are asked to complete a set of typical tasks while user experience researchers watch, listen, and take notes.



# Usability Test Methodology

# Participant Demographics

Macadamian recruited 20 practicing sonographers to participate in the usability test. Test participants were screened by third-party recruitment agencies for prior experience, brand attitude, and domain skills to mitigate bias.



100% specialize in abdominal sonography\*



Workload ranges:  
4 – 15 patients/ day  
20 – 50 patients/ week



Sonography experience:  
5 – 40 years



Hospital (12)  
Clinic (6)  
Both (2)




None of the participants were current Siemens ultrasound users

*\*The specialty of abdominal sonography was chosen in order to successfully recruit a large group of sonographers with similar demographics and experience in Spring 2015.*

# Usability Test Tasks

For the test, research participants were asked to perform a series of tasks that are commonly performed in an abdominal ultrasound exam. These were simulated by representative scanning tasks and were grouped into three categories.



## OPTIMIZE

E.g. Scan the aorta in color and pulsed wave Doppler modes



## ANNOTATE

E.g. Annotate three 2D images of the right kidney



## MEASURE

E.g. Measure the liver length and assign the measurement to a label





## Usability Test Space


The usability test was conducted at a third party user research facility in the Seattle, WA area equipped with a state-of-the-art usability testing lab.

Human models were involved in the study to act as patients during the usability tests so as to provide sonographers with real physiology to scan.

Both adult males and females, who were generally healthy, fairly easy-to-scan, and with no suspected abnormal anatomy, were recruited as models.

# Test Session Breakdown

~120 minutes



An introduction was provided first to participants which detailed the purpose of the study. Participants then watched a training video created by Macadamian that featured an experienced sonographer completing similar tasks with similar features, as included in the usability test tasks.

The usability test participants were instructed to watch the video first, followed by an opportunity to familiarize themselves with the system. Once training was complete, participants were given a break so that the usability test space could be set up for the test.

After completing all of the test tasks, a wrap-up session was conducted that included the System Usability Scale questionnaire.

# Usability Metrics (1)



The chosen metrics focused on effectiveness, efficiency and satisfaction, as these three elements define usability.

The international standard, ISO 9241-11, defines usability as:

The extent to which a product can be used by specified users to achieve specified goals with effectiveness, efficiency, and satisfaction in a specified context of use.

# Usability Metrics (2)

During the usability test, Macadamian's User Experience Researchers captured the data contained in the table below.

USABILITY ELEMENTS		CAPTURED DATA	
		QUANTITATIVE	QUALITATIVE
<b>EFFECTIVENESS</b>	Can users complete tasks and achieve goals with the system?	Success rates Number of errors	Think-aloud comments
<b>EFFICIENCY</b>	How much effort do users require to complete tasks?	Number of deviations	Think-aloud comments
<b>SATISFACTION</b>	What do users think about the system's ease of use?	Ease of use ratings	Think-aloud comments

Deviation Description(s)	Success	Number of Errors	Error Details	Comments	EOU Rating	Path Taken	Number of Deviations	Deviation Description(s)	Success	Number of Errors	Error Details	Comments	EOU Rating	Path Taken	Number of Deviations	Deviation Description(s)
-	1	0	-	It's easy, it's good	5	Scanned - Adjusted - Redo - Dual - Dual - Dual - Update - Freeze - It would be great if you can hit it again to switch onto the dual screen - freeze - update - gain	2.00	went in and out of dual, didn't set caliper properly at first	1	0	-	I like the ease of moving your calipers - just grab them and move them	4	ABC - SAG RIGHT KIDNEY - TRANS RIGHT KIDNEY ON tp - ABC - As soon as I hit ABC, it took my measurements off. My	0.00	-
-	1	0	-	Tp buttons are too close together	4	Right Kidney Sag on tp - Adjusted gain - Freeze - Update - Freeze - Freeze - ABC - Clear Screen - TRV - Clear Screen - Update - Gain	1.00	had to redo measurements	1	1	unfroze and had to do measurements again	Once you know a different machine, it would be pretty easy to use.	4	Right Kidney sag - freeze - dual - trv annotation	0.00	-
-	1	0	-	Its right there I didn't have to start over with the gb label.	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	5.00	had to redo image twice and accessed dual button	1	0	-	It would be nice if the image can take up more room on the monitor. Nice if you can measure off the image. The dual image was	3	ABC - Right Kidney Sag - Update - ABC - tb cines [lost a measurement] - ABC - length on tp - remeasured length - dragged	1.00	selected 'update'
-	1	0	-	3 Because I want to have post processing. 5 for annotations.	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	8.00	hit play by the track ball, freeze on and off (three times)	0	1	Hit update twice too fast	Expected when it came back to the original sector width. Clear everything. This doesn't make sense! Why do you have to	3	ABC - Right Kidney Sag on tp - moved trackball - TRV on tp	0.00	-
-	1	0	-	Nice to have a bigger keyboard picture on there	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	0.00	focus - zoom knob - that's interesting dual looks like it's multifunction - Freeze - Caliper - Clear Screen - Update - Gain	1	0	-	-	5	Refer to 4A	0.00	-
-	1	0	-	They are just right in front of me. Everything that I need.	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	1.00	selected dual and unfroze image	1	0	-	Just because when I unfroze I expected it to go to the next image - when I unfroze it it went back to single.	4	ABC - SAG RIGHT KIDNEY - MOVED Cursor - TRV right kidney - image	0.00	-
-	1	0	-	The ABC button is right there and then your labels come up.	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	2.00	had to redo image once in dual (twice)	1	0	-	It was easy to go into zoom mode - hard time figuring out how to zoom in on the dual screen. With the trackball it was left or	4	ABC - Tb to move - SAG RIGHT KIDNEY on tp - moved tb to other image - SAG - TRV RIGHT KIDNEY - image.	0.00	-
-	1	0	-	Umm I like the annotation is simple, it's just there. Supine and options.	4	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	2.00	had to redo image once in dual (twice)	1	0	-	Freezing the image and going to the right, and finding how to activate the other side. Sometimes it's just easier if it was	3	ABC - RIGHT KIDNEY SAG on tp - tb - TRV ON OTHER side	0.00	-
-	1	0	-	all I did was hit it and it automatically came up and I can see it.	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	4.00	had to redo image (four times)	0	1	Went to redo dual (i.e. same plane of the kidney) instead of	When you hit dual I couldn't get them both unless I hit dual live.	2	ABC - Tb - sag right kidney on tp - tb to other - TRV right kidney on tp image. [GB was still there]	0.00	-
Selected CBD annotation	1	2	Typed CBD twice	hitting the CBD below it and instead of deleting the last word it deleted the line. Oh it's there	3	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	1.00	had to redo image	1	0	-	After setting it I can move the arrow or caliper anywhere. It's easily allowing me to go through	5	Refer to 4A	0.00	-
-	1	0	-	Two clicks, turn on annotation, it's right there already.	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	1.00	went out of dual	1	0	-	Once I hit dual it wasn't intuitive how I got to the image	3	ABC - Delete - SAG Right kidney on tp - tb - TRV on tp - image.	0.00	-
-	1	0	-	this tp is really obvious. Only the stuff that I would need or use readily is on it.	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	5.00	had to redo image, went out of dual (twice), caliper (twice)	1	0	-	The freezing the image and saving and going to the next one is not convenient.	1	ABC - TRV using kb - tb - SAG using kb	0.00	-
-	1	0	-	You push the ABC button and the annotations come up and it's right there.	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	1.00	had to redo image	1	0	-	I couldn't figure out how to get the image on the right. Doesn't usually use Dual	4	Update twice - Kidney length on tp - ABC - RIGHT KIDNEY on tp - tb - SAG [on another line] - tb to move to the other image - TRV	0.00	-
-	1	0	-	I want - I want delete to do a clear screen function and take me back to home - set home at the bottom	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	1.00	redo measurements	1	1	Did set 2nd caliper and lost her measurement	Image problem in the machine... This limits the application of the dual function - I had to figure out what I am doing wrong - I had this	2	ABC - SAG RIGHT KIDNEY ON tp - MOVED Cursor - TRV right kidney - image	0.00	-
-	1	0	-	other than me having to push a button when it frozen. I couldn't find the set. Green cursor on the tp?	4	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	1.00	went out of dual	1	0	-	measurements: 2. I could have gotten the AP that would have taken up the screen - 15 images	2	ABC - RIGHT SAG KINEY on tp - tb - RIGHT TRV KIDNEY on tp	0.00	-
-	1	0	-	It was easy to find the annotation button and it popped up on the tp for me and it was easy to find on	5	Scanned - Update - Gain - Adjusted gain - Freeze - Update - Freeze - Update - Gain - Freeze - Update - Gain	5.00	had to redo image and measurements	1	2	measured on a single view and lost measurements.	4 for dual and for the measurement 2.	3	ABC - TRV using kb - tb - SAG using kb	0.00	-


# Usability Test Findings

# Effectiveness

task success, number of errors, qualitative feedback

## Average Task Completion Rate Across All Tasks

Error bars represent 95% confidence intervals



In terms of effectiveness, the results indicated that with the Siemens system **very few errors** were observed and tasks had **high completion rates**.




“I did like the annotations. I like how the ones that replace each other are all connected and there are no lines between them and the ones that are separated are the ones that don’t replace each other. I thought that’s cool because my current system doesn’t do that.”

*- Usability test participant 11*

# Efficiency

number of deviations, qualitative feedback



Regarding efficiency, **few deviations** were observed. For example, navigating to an incorrect touch panel screen or interacting unnecessarily with the interface.

A **large amount of positive qualitative feedback** was collected regarding efficiency of use.






“I really feel like an overarching concept for the design of this system was simplicity and elegance. The more simple it is, the easier it is to use. I did appreciate that there was a lot of streamlining”

*- Usability test participant 7*

# Satisfaction

ease of use ratings, qualitative feedback

The Siemens system scored **very high** for ease-of-use and also received a **large amount of positive qualitative feedback**.






"On my current machine, I use the keyboard a lot and I didn't even miss it."

*- Usability test participant 6*

# System Usability Scale

At the end of each usability test, the sonographers completed the industry standard System Usability Scale questionnaire to measure perceived system usability and learnability.



**For more information about this research and to learn more about how our test was conducted:**

**please email us at:** [info@macadamian.com](mailto:info@macadamian.com)

**please visit:** <http://www.macadamian.com/work-customer-experience/case-studies-software-application-developer/usability-testing-for-siemens-medical/>

m a c a d a m i a n

Thank you.