

# BeneFusion SP1/SP3/VP1/VP3 Vet

Veterinary Pump











inkodin E

Facebook



www.mindrayanimal.com service@mindrayanimal.com





## **BeneFusion SP3/VP3 Vet**

## SP3 Vet Syringe Pump

### Friendliness

- Compact and light, easy to carry and mount it
- One-hand operation and intuitive keypad
- Durable battery for emergency and patient transfer

### Safety

- High accuracy ≤±2% that is suitable for the use of small animals
- Fast-start function against start delay to ensure timely infusions
- Dynamic pressure display predicts occlusions in advance
- Anti-bolus function protects animal safety by releasing the pressure when an occlusion occurs

### Connectivity

- WiFi connection available
- Supports dock mounting

## VP3 Vet Infusion Pump

#### Friendliness

- Compact and light, easy to carry and mount it
- Durable battery for emergency and patient transfer
- Large screen with intuitive operation

### Safety

- Dynamic pressure display predicts occlusions in advance
- Anti-bolus function protects animal safety by releasing the pressure when an occlusion occurs
- Highly sensitive detectors:
- air sensors minimum 20µL
- pressure sensors minimum 75mmHg
- Sequential mode supports multiple infusions (piggyback)

### Connectivity

- WiFi connection available
- Supports dock mounting



# **BeneFusion SP1/VP1 Vet**

### SP1 Vet Syringe Pump

### Friendliness

- Compact and light, easy to carry and mount it
- One-hand operation and intuitive keypad
- Durable battery for emergency and patient transfer

### Safety

- High accuracy ≤±2% that is suitable for the use of small animals
- Fast-start function against start delay to ensure timely infusions
- Dynamic pressure display predicts occlusions in advance
- Anti-bolus function protects animal safety by releasing the pressure when an occlusion occurs

### Connectivity

- WiFi connection available
- Supports dock mounting

## VP1 Vet Infusion Pump

#### Friendliness

- Compact and light, easy to carry and mount it
- Durable battery for emergency and patient transfer
- Easy to operate and install IV sets

### Safety

- Dynamic pressure display predicts occlusions in advance
- Highly sensitive air bubble detection: minimum 50µL
- IP34 Protection

### Connectivity

- WiFi connection available
- Supports dock mounting



