


uMEC10 Vet

Quality Choice for Animal Primary Care



- 10.4-inch color LED backlight LCD
- Compact and portable design
- Vital parameter configuration
- Optimized workflow with friendly QuickKeys

Simplified cleaning experience

No animal hair stuck problem occurred with no fan inside, bring a simplified cleaning experience and long lifetime of monitor.



Accessories management

Unique cabinet design makes accessories management effective



Comprehensive Monitoring

With nearly 30 years of technical accumulation, uMEC10 Vet is the quality choice for animal primary care by providing accurate vital parameters measurements and animal-specific accessories.



Reliable measurements

NIBP: Specific algorithm improved with animal blood flow characteristics

SpO₂: Stable performance under weak perfusion or with interferences

ECG: Patented algorithm and Smart Lead ensure accurate monitoring

Specific Accessories



Crocodile clips with
comfort soft edge



Esophageal probe for
ECG & TEMP



Clip sensor for animal
tongue and ear



Multi-Size cuffs for
different animal sizes



LinkedIn



Facebook



Instagram



Search Mindray Animal Medical

www.mindrayanimal.com

service@mindrayanimal.com

P/N:ENG-uMEC10 Vet-210285X2P-20220913

©2022 Shenzhen Mindray Animal Medical Technology Co., LTD. All rights reserved.

